Brewers QCheck™ Kit Manual

WE BRING
IDEAS TO LIFE.

NATURAL INGREDIENTS
INGREDIENT SYSTEMS
INTEGRATED SOLUTIONS

Table of contents

	Page
1. Simple and reliable microbiological quality control	3
2. Contents: Brewers QCheck™ Kit	4
3. Storage	5
4. Intended purpose	5
5. Safety requirements	5
6. Quality control	6
7. Benefits	6
8. Materials needed for analysis – not included in kit	6
9. Important guidelines	8
10. Protocols	8
10.1. LMC: Water analysis	8
10.1.1. LMC sampling	8
10.1.2. LMC incubation	9
10.1.3. LMC results analysis	9
10.2. NBB®-B: Yeast and beer analysis	10
10.2.1. NBB®-B sampling	10
10.2.2. NBB®-B sample combination with culture medium	11
10.2.3. NBB®-B incubation	11
10.2.4. NBB®-B results analysis for yeast and beer	12
10.2.5. NBB®-B results analysis for the yeast and beer swabs	12
10.3. NBB®-B-Am: Production and dispensing unit hygiene (biofilm indicators)	13
10.3.1. NBB®-B-Am sampling for dispensing unit hygiene	13
10.3.2. Transferring the swab sample to the NBB®-B-Am tube	13
10.3.3. NBB®-B-Am sampling for production hygiene	13
10.3.4. Transferring the swab sample to the NBB®-B-Am tube	14
10.3.5. NBB®-B-Am incubation	14
10.3.6. NBB®-B-Am results analysis	15
11. Frequently Asked Questions (FAQs) on application and evaluating the results	16
11.1. LMC: Water analysis	16
11.2. NBB®-B: Yeast and beer samples	17
11.3. NBB®-B-Am: Production and dispensing unit hygiene	18
12. Appendix:	19
12.1. Brief overview: Brewers QCheck™ Kit – analyses	19
12.2. Glossary	20
13. Supply sources	21
14. Information	22

1. Simple and reliable microbiological quality control

One of the main tasks in breweries is conducting microbiological beer fermentation in a clean and hygienic way. However, beer spoiling microorganisms can inhibit fermentation or have a negative impact on the beer's flavours. Therefore, microbiological control of the beer brewing processes is extremely important and critical.

The Brewers QCheck™ Kit makes it possible to perform microbiological control on all samples of the brewing process in a simple, fast and comprehensive way. Different raw materials used for the brewing process, such as water and yeast, can be analysed for microbiological contamination before fermentation. Beer spoiling bacteria can be detected in the end product, i.e. the beer, and the hygiene of the surrounding area, for example within the production site and dispensing unit, can be accurately monitored.

All microbiological detection is easily and reliably visible through a colour change, even without any previous knowledge. Consequently, the Brewers QCheck™ Kit allows the entire beer brewing process to be microbiologically controlled, thus providing the perfect introduction to microbiological quality control.

2. Content: Brewers QCheck™ Kit

The Brewers QCheck™ Kit, Article No. 2.11753.244, contains the following individual components from the Döhler Microsafety Design® (DMD®) portfolio which are applied to the microbiological analysis of raw materials used during the brewing process, of the final beer and hygiene of production and dispensing units.

Döhler Microsafety Design® Product	LMC	NBB®-B Tubes	NBB®-B-Am Tubes	Swabs, sterile**
Description	For water analyses in accordance with DEV*: Ready-to-use lactose peptone concentrate for the detection of faecal indicator bacteria	For yeast analyses of propagation, pitching, assimilation and harvesting yeast, as well as beer: Ready-to-use NBB®-B Broth for detecting beer spoiling microorganisms by a colour change	For hygiene monitoring of production and dispensing units: Ready-to-use NBB®-B-Am Broth for the detection of biofilm indicator and obligate beer spoiling microorganisms by a colour change	For sample collection of yeast beer and beer containing yeast from production, contact surfaces and dispensing units.
Target germ(s)	<i>E. coli</i> /Coliform bacteria	Beer spoiling microorganisms: Lactobacilli, Pediococci, Pectinatus and Megasphaera	Biofilm indicator microorganisms and beer spoiling microorganisms, e.g.: Acetic acid bacteria, Lactic acid bacteria, Pectinatus and Megasphaera	–
Number/kit	2 bottles with 50 ml of LMC concentrate in a carrier bag	1 polystyrene box with 20 tubes of 10 ml NBB®-B	1 polystyrene box with 20 tubes of 10 ml NBB®-B-Am	3 x 10 sterile swabs (individually packaged)

*DEV = German standard methods for drinking water analysis

** Also available without the kit with 100 pc. per pack (see 13. Supply sources; p. 21)

3. Storage

The kit and its individual components must be stored in a dry place away from the light at temperatures ranging between 4 - 8°C. The sterile swabs must be stored in a dry place at room temperature. The products shall not be frozen. An unopened kit can be used in the conditions outlined above until the best-before date expires. The expiration date can be found on the label on the packaging for the Brewers QCheck™ Kit. It may well be that the expiration dates of the individual components exceed the expiration date specified on the Brewers QCheck™ Kit (see product label) and can be used for a longer period of time.

When you dispose of the packaging, please note the batch number on the Brewers QCheck™ Kit (see label) for traceability.

4. Intended purpose

The purpose of the Brewers QCheck™ Kit is to analyse contaminations found in the beer brewing process as well as the hygiene of production and dispensing units. When using the kit, we recommend that you pay attention to the safety requirements outlined below and work particularly cautiously and hygienically to prevent secondary contamination during sampling.

5. Safety requirements

The basic rules for microbiological operations shall be observed when using this product. This includes using clean equipment such as lab coats, safety goggles and gloves. Such equipment contributes to your own personal safety and prevents contamination by the user. The Material Safety Data Sheets (MSDS) for the individual kit components can be found at www.doehler-dmd.com. No component of the Brewers QCheck™ Kit is suitable for consumption.

Special care is required when using Bunsen burners while wearing latex gloves. Always maintain a sufficiently safe distance from the flame. If the gloves catch fire, it can result in extremely severe injuries.

Fire hazard:

Never use alcohol for sterilisation or decontamination purposes while working with an open fire (Bunsen burner).

6. Quality control

The quality of the Brewers QCheck™ Kit and its individual components is carefully examined. In particular, the culture media's functionality is specifically tested, proved and confirmed using test microorganisms.

7. Benefits

The Brewers QCheck™ Kit contains ready-to-use culture media, eliminating the need for extensive and costly preparations such as weighing and mixing the media components, adjusting pH values and subsequent autoclaving. This saves you both time and money, thus allowing you to immediately focus your efforts on carrying out analyses of your yeast, beer and brewing process.

8. Materials needed for analysis – not included in kit

The aim of this product is to be able to conduct microbiological quality control checks with just a minimum of laboratory equipment. However, the following items are crucial for conducting the experiments:

Labware	Use
<p>Spray bottle containing alcohol solution used for sterilising surfaces/disposable gloves.</p> <p>(e.g. 70% ethanol/ 30% distilled water)</p>	<p>Spray on work surfaces and leave for at least 5 minutes, then wipe the surfaces with a paper towel.</p> <p>Before carrying out microbiological operations, decontaminate disposable gloves with a quick spray.</p> <p></p> <p>Warning: Flammable! Never use in the presence of an open flame!</p>

<p>Portable Bunsen burner</p>	<p>Due to their flames, Bunsen burners create a sterile work area (approx. 50 cm around the flame) in which microbiological experiments can be carried out. The heat kills bacteria in the environment. A portable Bunsen burner is also a useful tool for sterile sampling.</p>
<p>A sterile, serological pipette (2 ml) with pipetting device (e.g. Peleus ball or plastic pipetting device)</p>	<p>Serological pipettes can be purchased sterile. They are suitable for transferring sterile liquids. A large number of samples can be extracted using just the swab included. However, it is essential to use the pipette for clear beer samples.</p> <p style="text-align: center;"></p> <p>Warning: never pipette fluids with your mouth! Use a pipetting aid at all times!</p>
<p>Test tube racks</p>	<p>Rack for NBB®-B and NBB®-B-Am Tubes</p>
<p>Incubator with thermostat</p>	<p>For incubating samples at certain temperatures</p> <p>Useful tip: Some rooms within the brewery may have stable temperature conditions. These rooms are equipped with a thermometer for temperature measurements and can serve as a practical alternative to an incubator</p>
<p>Disposable gloves</p>	<p>Disposable gloves prevent secondary contamination by the user</p>
<p>Sterile 50 ml sample tubes for beer sampling</p>	<p>For example Döhler TransFast® Tubes, Article no. 2.04730.001. Equivalent products by other manufacturers can also be used.</p>

9. Important guidelines

Clean microbiological working is the basis of all analyses. It is highly recommended to follow the work steps described (see 10. Protocols, p. 8) in order to obtain reliable results.

No culture media are suitable for consumption.

Following analysis, the incubated culture media should be deactivated before being disposed. Autoclaving the culture media or disposing of them as hazardous waste is recommended.

10. Protocols

10.1. LMC: Water analysis

10.1.1. LMC sampling

Briefly flame the water tap if the system is made of metal. For plastic components, use the alcohol solution for decontamination.

Let the water flow from the tap for approx. 20 - 30 seconds.

Open the LMC bottle.

Remove the lid from the bottle neck just enough to be able to fill up the bottle.

Fill the LMC bottle with approx. 250 ml of water (for filling levels, please refer to 10.1.3 results analysis). Please do not fill the bottle up completely!

Pass the bottle neck through the flame using the portable Bunsen burner and seal the filled bottle.

10.1.2. LMC incubation

Incubate the sealed LMC bottle at 37 ± 2 °C in an incubator for 2 days.

10.1.3. LMC results analysis

Once incubation is complete, the colour indicator of the medium will turn yellow if *E.coli* or Coliform microorganisms are present. The solution will also appear cloudy and gas formation can be observed.

Result: The water is contaminated and should not be used in the brewing process. Repeat the experiment to confirm the results.

10.2. NBB[®]-B: Yeast and beer analysis

For yeast and cloudy beer, it is possible to use the swabs provided and a sterile serological pipette (not available with the kit) for sampling. For clear beer samples, please use the serological pipette.

10.2.1. NBB[®]-B sampling

Yeast:

Remove 0.5 - 1 ml of yeast from the yeast culture using a sterile serological pipette.

Alternative: It is possible to extract thick, viscous yeast using the sterile swab provided.

Beer:

Remove a beer sample (approx. 50 ml, see Useful tip: beer samples) under sterile conditions and transfer this to a sterile sample flask.

A) Sampling with a pipette

Remove 0.5 - 1 ml of (green) beer in the sterile sample flask with the beer sample (50 ml) or alternatively from a bottle of beer, using a sterile serological pipette.

B) Sampling with a swab:

Leave the beer sample (50 ml) or beer bottle for 24 hours at 4°C to enable the formation of sediment at the bottom. Decant the liquid. Remove a sample from the sediment using a sterile swab.

Useful tip: beer samples

The majority of samples are taken from a tank using a sampling tap, or are alternatively taken from a barrel using a fitting. When carrying out this procedure, it is important that the sampling tap is decontaminated by using alcohol solution or a Bunsen burner. Valves/ fittings can be immersed in a suitable disinfectant beforehand. As it is not uncommon for severe contaminations to be found in the areas of the sampling tap, it is important to leave the beer running for a short while (approx. 1 second). Then take a sample.

10.2.2. NBB[®]-B sample combination with culture medium

Yeast:

Transfer the yeast sample into an NBB[®]-B Tube.

Pass the tube through the flame and seal it.

Alternative: Insert the swab into the NBB[®]-B Tube.

Break off the upper part of the wooden shaft, flame the tube and seal it.

Beer:

A) Sampling with a pipette

Fill an NBB[®]-B Tube with 0.5 - 1 ml of the beer sample.

Flame the NBB[®]-B Tube and seal it.

B) Sampling with a swab:

Insert the sterile swab with the yeast sediment into an NBB[®]-B Tube.

Break off the upper part of the wooden shaft and seal the tube.


10.2.3. NBB[®]-B incubation

Yeast:

Incubate the sealed NBB[®]-B Tube at $28 \pm 2^{\circ}\text{C}$ for 5 days.

Beer:

Incubate the sealed NBB[®]-B Tube at $28 \pm 2^{\circ}\text{C}$ for 5 days.

10.2.4. NBB®-B results analysis for yeast and beer

Once incubation is complete, the NBB®-B Tube colour will change from red to yellow if beverage spoiling microorganisms have been positively detected. The solution will also appear cloudy and gas formation may be observed in some cases.

Result: The yeast is contaminated with beer spoiling microorganisms and should not be used in the brewing process.

The beer is contaminated with beer spoiling microorganisms and should not be bottled or delivered.

If possible, repeat the experiments to confirm the results. If the beer is contaminated, it is necessary to carry out further analyses in order to determine its potential risk.

10.2.5. NBB®-B results analysis for the yeast / beer swabs

Once incubation is complete, the NBB®-B Tube colour will change from red to yellow if beverage spoiling microorganisms have been positively detected. The solution will also appear cloudy and gas formation may be observed in some cases.

Result: The yeast is contaminated with beer spoiling microorganisms and should not be used in the brewing process.

The beer is contaminated with beer spoiling microorganisms and should not be bottled or delivered.

If possible, repeat the experiments once again to confirm the result. If the beer is contaminated, it is necessary to carry out further analyses in order to determine its potential risk.

10.3. NBB®-B-Am: Production and dispensing unit hygiene (biofilm indicators)

10.3.1. NBB®-B-Am sampling for dispensing unit hygiene

Leave the beer dispensing unit to run for approx. 1 - 2 seconds.

Hold the sterile swab in the stream of beer for approx. 15 seconds.

10.3.2. Transferring swab sample to the NBB®-B-Am Tube

Open the NBB®-B-Am Tube only briefly to insert the swab. Carefully break off the wooden shaft in the test tube.

Flame the NBB®-B-Am Tube and seal it.

10.3.3. NBB®-B-Am sampling for production hygiene

Specify approx. 7 - 12 defined control points within the filling area where you carry out routine sampling (e.g. valves, star connectors, casings, stamps, organs (guides), etc.). They should be monitored on a weekly basis.

Take a test tube rack and place a sufficient number of NBB®-B-Am Tubes inside. You need an NBB®-B-Am Tube and a sterile swab for every control point.

Take them to the control point in the production.

If the sample point is already moistened, you can immediately take a sample using the swab.

If the sample point is not moist, open the NBB[®]-B-Am Tube quickly and moisten the swab using the media solution. Take a sample of the control point using the moistened swab.

Warning: Never moisten the swab again if it has been taken from the tube; this may cause a secondary contamination in the tube.

10.3.4. Transferring swab sample to the NBB[®]-B-Am Tube

Open the NBB[®]-B-Am Tube only briefly to insert the swab.

Carefully break off the wooden shaft in the test tube. Flame the NBB[®]-B-Am Tube and seal it.

Repeat this procedure for each control point.


10.3.5. NBB[®]-B-Am incubation

Incubate all NBB[®]-B-Am Tubes with a swab in an incubator at $27 \pm 2^{\circ}\text{C}$ for 3 days maximum.

Warning: Please follow the incubation timeframe; contaminations that occur after 3 days are not linked to biofilm indicator microorganisms.

10.3.6. NBB®-B-Am results analysis

The sample points are analysed after 3 days' incubation.


Strong biofilm formation occurs if the NBB®-B-Am Tube turns yellow after just one day and becomes very cloudy. Gas formation may also be observed in some cases.

All samples are analysed to observe the overall hygiene condition.

Results: The hygiene condition is “perfect” for less than 20 % of positive findings. The test results of an interval must be documented.

If the number of positive findings continuously increases, i.e. more than 20 % of positive samples are found for two to three sample intervals (see example in figure below), the production facility and devices need to be cleaned. The risk of beer contamination can be prevented.

Amount of positive NBB®-B-Am Tubes

Biofilm development in a production facility using the positive findings of NBB®-B-Am Tubes (blue). The 20% safety threshold is presented as a black broken line. Production hygiene under this threshold is “perfect”. If the safety threshold has been passed, biofilm development is shown after cleaning (green) and no cleaning (red). If you identify values above the 20% safety threshold again, microorganisms may be enriched in biofilms, thus causing secondary or scattered infections.

11. Frequently asked questions (FAQs) on application and evaluating the results

11.1. LMC: Water analysis

Bottle running over

Cause: You have poured too much sample into the bottle. The sample volume is 250 ml.

Solution: Please check the sample volume and do not incubate the overfilled bottle.

Is it possible to differentiate between *E. coli*/Coliform microorganisms?

Cause: Qualitative detection of *E. coli*/Coliform microorganisms is carried out using the product provided. It is not possible to identify the microorganisms.

Solution: Other culture media can be used for identifying the microorganisms.

No bacterial growth at 28°C.

Cause: *E. coli*/Coliform microorganisms are incubated at $37 \pm 2^\circ\text{C}$.

Solution: Incubate the samples in LMC at $37 \pm 2^\circ\text{C}$.

Turbidity with no change in colour

Cause: Contaminations not with *E. coli*/Coliform microorganisms.

Solution: If necessary, conduct a subsequent analysis on the contamination using other culture media for identification.

11.2. NBB[®]-B: Yeast and beer samples

Turbidity with no change in colour

Cause: If the culture medium turns cloudy after 5 days incubation without colour change, this indicates contamination with no beer spoiling microorganisms.

Solution: Please repeat the sampling and follow the outlined work steps. Ensure that microbiological work is conducted hygienically.

Identifying beer spoiling microorganisms

Cause: It is not possible to identify beverage spoiling microorganisms using the Brewers QCheck[™] Kit.

Solution: Further analysis techniques can be applied for the identification of microorganisms.

Yeast in tube

Cause: Shortly after adding the yeast in larger quantities, temporary growth may result in a colour shift in the NBB[®]-B Tube. However, the yeast's growth will be inhibited after a short period. The yeast may also cause turbidity in the tube.

Solution: After 5 days, only beer spoiling microorganisms are detected in the analysis. Please follow the exact analysis times. If the solution turns cloudy, please check for a change in colour as this indicates the presence of beer spoiling microorganisms.

A colour change immediately after adding beer

Cause: You have poured too much beer into the NBB[®]-B Tube. This will immediately change the colour indicator.

Solution: Please follow the outlined process steps. Do not add more than 0.5 - 1 ml beer per tube with the pipette.

Turbidity immediately after adding beer

Cause: The tube may appear slightly cloudy after adding the cloudy beer sample or green beer.

Lösung: The analysis of the tube is based on the colour change. There may also be gas formation.

11.3. NBB®-B-Am: Production and dispensing unit hygiene

Turbidity with no change in colour after 3 days of incubation

Cause: Turbidity may appear after 3 days of incubation. The microorganisms detected are not linked to biofilm formation and are in most cases harmless and ubiquitous.

Solution: Please follow the incubation times and carry out the analysis after 3 days.

| 12. Appendix

12.1. Brief overview: Brewers QCheck™ Kit – analyses

DMD® Product/ sample	Sample volume	Incubation conditions	Analysis
LMC: Water analysis	250 ml of water in LMC bottle	Temp.: 37 ± 2°C Time: 2 days	Positive: Colour change, purple to yellow, turbidity, (gas formation) Negative: No colour change, no turbidity
NBB®-B Tube: Yeast	0.5 - 1 ml of yeast or 1 swab sample	Temp.: 27 ± 2°C Time: 5 days	Positive: Colour change, red to yellow, clear turbidity, (gas formation) Negative: No colour change, no additional turbidity
NBB®-B Tube: Clear beer	0.5 - 1 ml of beer sample	Temp.: 27 ± 2°C Time: 5 days	Positive: Colour change red to yellow, turbidity, (gas formation) Negative: No colour change, no turbidity
NBB® B Tube: (Yeast) cloudy beer	0.5 - 1 ml of the cloudy beer sample or swab with sediment	Temp.: 27 ± 2°C Time: 5 days	Positive: Colour change red to yellow, (gas formation) Negative: No colour change, no additional turbidity
NBB®-B-Am Tube: Dispensing units hygiene	1 swab sample in stream of beer for 15 seconds	Temp.: 27 ± 2°C Time: 3 days maximum	Positive: Colour change red to yellow, turbidity, (gas formation) Negative: No colour change, no turbidity
NBB®-B-Am Tube: Production hygiene	7 - 12 swab samples	Temp.: 27 ± 2°C Time: 3 days maximum	Positive: Colour change red to yellow,turbidity, (gas formation) Negative: No colour change, no turbidity

12.2. Glossary

Identification	Identification is the differentiation of bacteria on the basis of their genus and species names; e.g. <i>Escherichia coli</i> . Identification is the determination of microorganisms.
Incubation	The process of incubating a culture media under defined conditions and parameters (e.g. in terms of time and temperature)
Peleus ball	A rubber ball with valves used for sucking in and discharging liquids by means of glass or plastic pipettes.
Qualitative detection	Detection based on a yes/no statement; this means that microorganisms are or are not contained in the sample.
Quantitative detection	Detection based on a yes/no statement; this means that microorganisms are or are not contained in the sample, the result can be counted.
Serological pipette	Pipette usually made out of plastic or glass used for transferring liquids. The scale featured on the pipette displays the volume. Plastic pipettes are usually available in a sterile condition and are individually packaged.
Sterile	Sterile means free from bacteria or living microorganisms.

13. Supply sources

The company Döhler provides you with different solutions for quality control checks:

Article no.	Product	Packaging unit
2.11753.244	Brewers QCheck™ Kit	1 kit

The products contained in the kit are also available in individual larger packaging units:

Article no.	Product	Packaging unit
2.04713.700	LMC	9 x 50 ml in bottle (cardboard box)
2.04723.646	NBB®-B Tubes	20 x 10 ml in tube (box)
2.04706.646	NBB®-B-Am Tubes	20 x 10 ml in tube (box)
2.04725.244	Sterile swab (sterile)	100 pc./ package

Other Döhler Microsafety Design® products:

Article no.	Product	Packaging unit
2.04730.001	TransFast® Gel Tube (sterile)	100 pc./ package

All the labware that is not available in the kit can be obtained from a laboratory supplier.

| 14. Information

Trademarks

NBB®, Döhler Microsafety Design® und Brewers QCheck™ are all trademarks by Döhler GmbH that are registered and protected worldwide.

Further reading

Back, W.: Colour Atlas and Handbook of Beverage Biology, Fachverlag Hans Carl, Nuremberg, 2005.

Back, W. (publisher): Selected chapters of brewing technology, second updated edition, pg. 319 seq., Fachverlag Hans Carl, Nuremberg, 2008.

Distribution

You can find more information about our Döhler Microsafety Design® product portfolio on the website www.doehler-dmd.com. Email contact: dmd@doehler.com

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